


# NIPPY S+



## Bi-Level positive pressure ventilator

for use with Nasal Mask or Nasal Pillows

### Optional INTERNAL BATTERY

*Run times up to 12 hours*

The **Nippy S+** is a pressure controlled positive pressure ventilator employing state of the art electronics, to give consistent and accurate support for patients requiring non invasive ventilation via a nasal mask or nasal pillows.

The large colour LCD screen gives unrivalled clarity when setting up or reading stored data.

Intuitive software allows for easy setting up, with on-screen help available for each function.

### Pressure Support Ventilation

Settings and airway measurements are displayed on-screen with more advanced features accessed via an easy to follow menu.

Compliance data may be viewed on-screen.

Full event log may be downloaded to a PC.

The **Nippy S+** is equipped with a universal mains power inlet, a dc power inlet and an

Tel: 01789 293460 Sales.  
Fax: 01789 262470


01789 721577 Technical Support  
[www.nippyventilator.com](http://www.nippyventilator.com)

internal battery. Weighing only 3.6kg,  
world-wide travel has never been easier.

Tel: 01789 293460 Sales.  
Fax: 01789 262470

01789 721577 Technical Support  
[www.nippyventilator.com](http://www.nippyventilator.com)

# Nippy S+



The Nippy S+ is intended for use by patients who require respiratory support via nasal mask or nasal pillows.

For the treatment of:-  
Chronic Obstructive Pulmonary Disease (COPD)

Decreased drive to breathe (Central Hypoventilation)

The Nippy S+ is equally suited to hospital or domicilliary use.

## Technical Specifications

<b>Supply Voltage</b>	100 - 240 Vac
<b>Dimensions cm (LxWxH)</b>	30 x 22 x 13
<b>Weight</b>	3.6 kg
<b>With battery</b>	4.5kg
<b>Pressure Range IPAP</b>	3 – 38 cm H <sub>2</sub> O
<b>Pressure Range EPAP</b>	3 – 20 cm H <sub>2</sub> O
<b>Fixed Back-up Rate</b>	12BPM
<b>Ti</b>	1.5 sec
<b>Optional Internal battery</b>	24Vdc 6.8 Ahr
<b>Running time</b>	4 -12 hours
<b>Optional External battery</b>	24Vdc 6.8 Ahr
<b>Running time</b>	4 -12 hours
<b>Alarms</b>	Mask off, Power Failure

## Features

- ◆ Pressure Support – Spontaneous breathing augmented by the ventilator
- ◆ User friendly intuitive software
- ◆ Large colour LCD display, clearly shows all settings.
- ◆ Lockable settings.
- ◆ Alarms with easily understood on-screen messages.
- ◆ On-screen compliance data.
- ◆ Full event log showing all adjustments and alarm events may be downloaded to a PC.
- ◆ Universal mains input and internal battery operates anywhere in the world.
- ◆ Adjustable fast response flow triggers.
- ◆ Can be connected to Nurse call systems using the optional nurse call adapter.

Tel: 01789 293460 Sales.  
Fax: 01789 262470

01789 721577 Technical Support  
[www.nippyventilator.com](http://www.nippyventilator.com)